




MB GB-2202

Single piece molded Glovebox

INERTGAS TECHNOLOGY



- Unique single piece molded design of box body
- No glue or clamps on the body of the box
- Easy to open or remove clamped front window
- Variety of glove options, standard is Neoprene
- Ball valves for gas purge and exhausting of the box
- Side door \varnothing 290 mm
- Bench top system

MB GB-2202

Standard Glovebox made of acrylic glass **MB GB-2202** with removable front window and side door $\varnothing 290$ mm

Article no. 7050148

Basic equipment	Technical data
<p>Glovebox</p> <ul style="list-style-type: none">• Single piece molded housing• Removable front window made of acrylic glass, with 2 glove ports• Pair of gloves made of Neoprene, size 8, lh/rh• Rubber ring fasteners for securing of gloves• Window clamps• Ball valve G 1/4", with 2 hose connections, for gas inlet and outlet of the Glovebox <p>Basics</p> <ul style="list-style-type: none">• The Glovebox can be operated either:<ul style="list-style-type: none">- at a negative pressure down to: -2,5 mbar- at a positive pressure up to: 2,5 mbar	<p>Glovebox properties</p> <ul style="list-style-type: none">• Length: 950mm• Depth below: 540mm• Depth above: 420mm• Height: 650mm • Volume of Glovebox: 300l • Weight: 21kg
<p>Side door:</p> <ul style="list-style-type: none">• Horizontal opening door $\varnothing 290$mm, made of aluminium	<p>Options:</p> <ul style="list-style-type: none">• Electrical power feedthrough, single and multiple internal connections• Low voltage feedthrough• Stainless steel floor pan• Internal and external covers for glove ports• Automated pressure regulation• Lighting options• Antechamber



M. Braun Inertgas-Systeme GmbH Head-Office Germany • Dieselstr. 31, D-85748 Garching • Phone: +49(0) 89-3 26 69-0
Fax: +49(0) 89-3 26 69-105 • E-Mail: info@mbraun.de

M. Braun Inc. USA • 14 Marin Way, Stratham, NH 03885 • Phone: +1-603-773-9333 • Fax: +1-603-773-0008 • E-Mail: info@mbraunusa.com

M. Braun Inertgas Systems (Shanghai) Co., Ltd. China • 828 Xin Jinqiao Road • Pudong Shanghai 201206, P.R.C.
Phone: +86-21-50320257 • Fax +86-21-50320229 • E-Mail: info@mbraunchina.com

M. Braun Service Japan • Yoshikawa building, 3rd floor • 6-32 Honcho, Wako-shi • Saitama 351-0114 • Phone: +81-48-452-6367 • Fax: +81-48-452-6368

Internet www.mbraun.de • www.mbraunusa.com • www.mbraunchina.com